Cochlear™

# Baha® 5 SuperPower Sound Processor

## User manual part A

English Svenska Dansk Suomi Norsk **ZONE 3**

3 15 27 39 51  
GB SV DA FI NO

*Hear now. And always*

## Sound processor type designations for models included in this User Manual are:

FCC ID: QZ3BAHA5SUP, IC: 8039C-BAHA5SUP, IC model: Baha®5 SUP.

### Statement:

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from the one in which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- Changes or modifications can void the user's authority to operate the equipment.

### Intended use

The Cochlear™ Baha® 5 SuperPower uses bone conduction to transmit sounds to the cochlea (inner ear). It is indicated for people with conductive hearing loss, mixed hearing loss and single sided sensorineural deafness (SSD). Furthermore it is indicated for bilateral and paediatric recipients. Fitting range up to 65 dB SNHL. It works by combining a sound processor and a small titanium implant that is placed in the skull behind the ear. The skull bone integrates with the titanium implant through a process called osseointegration. This allows sound to be conducted via the skull bone directly to the cochlea, which improves hearing performance. The sound processor can be used together with the Baha Softband.

The fitting is to be done either at a hospital, by an audiologist, or in some countries, by a hearing care professional.

### List of countries:





Not all products are available in all markets. Product availability is subject to regulatory approval in the respective markets.

The products are in compliance with the following regulatory requirements:

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Equipment includes RF transmitter


13

14


15

16


17




18

19

20

# Welcome

Congratulations on your choice of the Cochlear™ Baha® 5 SuperPower Sound Processor. You are now ready to use Cochlear's highly advanced bone conduction sound processor, featuring sophisticated signal processing and wireless technology.

This manual is full of tips and advice on how to best use and care for your Baha sound processor. By reading this manual and then keeping it handy for future reference, you will ensure that you get the most benefit out of your Baha sound processor.

## Key to device *See figure* **1**

1. Sound Processing unit
2. Actuator unit
3. Cable
4. Microphones
5. Battery
6. Ear hook
7. Plastic snap connector
8. Buttons
9. Visual indicator

**Note on figures:** The figures included on the cover correspond to information specific to this model of sound processor. Please reference the appropriate figure when reading. The images shown are not to scale.

## Introduction

To ensure optimal performance, your hearing care professional will fit the processor to suit your needs. Be sure to discuss any questions or concerns that you may have regarding your hearing or use of this system with your hearing care professional.

## Warranty

**The warranty does not cover defects or damage arising from, associated with, or related to the use of this product with any non-Cochlear processing unit and/or any non-Cochlear implant. See the "Cochlear Baha Global Limited Warranty card" and the "Addendum to Warranty Card" for more details.**

## **Contacting Customer Service**

We strive to provide you with the best possible products and services. Your views and experiences with our products and services are important to us. If you have any comments you would like to share, please contact us.

### **Customer Service – Cochlear Americas**

13059 E. Peakview Avenue, Centennial  
CO 80111, USA

Toll free (North America) 1800 523 5798

Tel: +1 303 790 9010, Fax: +1 303 792 9025  
E-mail: [customer@cochlear.com](mailto:customer@cochlear.com)

### **Customer Service – Cochlear Europe**

6 Dashwood Lang Road, Bourne Business Park,  
Aldershot, Surrey KT15 2HJ, United Kingdom

Tel: +44 1932 26 3400, Fax: +44 1932 26 3426  
E-mail: [info@cochlear.co.uk](mailto:info@cochlear.co.uk)

### **Customer Service – Cochlear Asia Pacific**

1 University Avenue, Macquarie University,  
NSW 2109, Australia

Toll free (Australia) 1800 620 929  
Toll free (New Zealand) 0800 444 819

Tel: +61 2 9428 6555, Fax: +61 2 9428 6352  
or Toll free Fax 1800 005 215

E-mail: [customerservice@cochlear.com.au](mailto:customerservice@cochlear.com.au)

## Key to symbols

The following symbols will be used throughout this document.

Please refer to the list below for explanations:

	"Consult instructions for use" or "Operating instructions"		Bluetooth® Smart
	"Caution" or "Caution, consult accompanying documents"		Manufacturer
	Audible signal		Catalogue number
	CE-mark		Date of manufacture
	Interference risk		By prescription
	Batch Code		Temperature limitation
	Serial Number		Humidity limitation
	Refer to instructions/booklet Note: Symbol is blue		Recyclable material
	ACMA (Australian Communications and Media Authority) compliant		Separate disposal of electronic device
	Made for iPod, iPhone, iPad		Do not re-use
			
		Product name, CE mark and ACMA symbol	

# Using your sound processor

The buttons on your sound processor let you adjust your volume level, choose from your pre-set programs, and enable/disable wireless streaming. You can choose to enable the audible and visual indicators to alert you of changes to the settings and processor status.

Your sound processor is programmed to be used as either a left or right-sided device. Your hearing care professional will have marked your processor(s) with a L or R indicator. If you are a bilateral user, changes you make to one device will automatically apply to the second device.

## On/off *See figure 2*

- 1 2 Turn on your sound processor by connecting the battery module.
- 3 Turn off your sound processor by disconnecting the battery.

When your sound processor is turned off and then back on again, it will return to the default setting (program one).

## Status indicator *See figure 3*

Your sound processor is equipped with audible and visual indicators. For an overview of the audible and visual indicators, refer to the table at the back of this section.

**Your hearing care professional can disable the audible and/or visual indicators if you prefer.**

## Change program *See figure 4*

Together with your hearing care professional you will have selected up to four pre-set programs for your sound processor:

- Program 1: \_\_\_\_\_
- Program 2: \_\_\_\_\_
- Program 3: \_\_\_\_\_
- Program 4: \_\_\_\_\_

These programs are suitable for different listening environments. Ask your hearing care professional to fill in your specific programs.

To switch programs, press and hold the upper button on your sound processor. If enabled, audible and visual indicators will let you know which program you are using:

- Program 1:** 1 beep and 1 flash
- Program 2:** 2 beeps and 2 flashes
- Program 3:** 3 beeps and 3 flashes
- Program 4:** 4 beeps and 4 flashes

## Adjust volume

*See figure 5*

- 1 Press and release the upper button to increase volume.
- 2 Press and release the lower button to decrease volume.

**You can adjust the volume level with the optional Cochlear Baha Remote Control, the Cochlear Wireless Phone Clip or an iPhone, iPad or iPod touch (see Made for iPhone section).**

## Wireless accessories

You can use Cochlear Wireless accessories to enhance your listening experience. Ask your hearing care professional to learn more about your options or visit [www.cochlear.com](http://www.cochlear.com).

**To activate wireless audio streaming,** press and hold the lower button until you hear a melody. *See figure 6*

**To end wireless audio streaming,** press and hold the upper button. The sound processor will return to the previous program.

## Flight mode *See figure 8*

When boarding a flight, wireless functionality must be deactivated because radio signals may not be transmitted during flights. To disable wireless operation:

- 1 Turn off the sound processor by disconnecting the battery.
- 2 Press and hold the lower button and connect the battery at the same time.

To deactivate flight mode, turn the sound processor off and back on again by disconnecting and connecting the battery.

## Made for iPhone MFi

Your sound processor is a Made for iPhone (MFi) hearing device. This allows you to control your sound processor and stream audio directly from iPhone, iPad or iPod touch. For full compatibility details and more information visit [www.cochlear.com](http://www.cochlear.com).

## Talk on the phone *See figure 7*

For best results, use a Cochlear Wireless Phone Clip or stream the conversation directly from your iPhone. When using a normal handheld phone, place the receiver near the microphone inlet of your sound processor instead of near your ear. Ensure the receiver does not touch the sound processor as this may result in feedback.

## Batteries

There are two different batteries for the Baha 5 SuperPower: the Standard Rechargeable Battery Module and the Compact Rechargeable Battery Module. Ask your hearing care professional for more information about the choice of battery.

## Battery life

Batteries should be replaced as needed just as you would with any other electronic device. Battery life varies according to the programs used each day, the sound environment use of wireless accessories, and the size of battery.

The rechargeable battery's lifespan is at least 365 charge cycles. A completely empty rechargeable battery will take approximately four hours to recharge.

Rechargeable batteries may take longer to fully recharge with age. To get the longest life from the rechargeable batteries, always recharge before use.

When there is approximately one hour of battery power remaining, you will hear an audio indicator comprising of two sets of four short beeps, and the visual indicator will flash repeatedly. At this time you may experience low amplification. If the battery runs down completely, the sound processor will beep and then stop working.

## Replace the battery

### Remove the battery *See figure 2*

Twist the battery module to remove it from the processing unit.

### Attach the battery *See figure 2*

Align the battery module to the processor socket to fit the parts together. Twist the battery module as shown to attach the parts. The processor will turn on automatically.

### Lock/Unlock the battery (tamper-proof)

*See figure 9*

- 1 To increase tamper resistance, push the lock to the far left with your battery cover locking tool to lock the battery module.
- 2 Push the lock to the far right to unlock the battery module.

### Charge battery modules *See figure 10*

- 1 Angle and fit the rechargeable battery module into a socket on the Cochlear Nucleus® Battery Charger.
- 2 Twist clockwise to connect.
- 3 Connect the battery charger to the power adaptor and plug into mains power.

## Care for your battery charger

### Every day

Check that your battery charger is clean. If you notice any dust or dirt:

- 1 Unplug the power adaptor and remove any battery modules.
- 2 Hold the battery charger upside down and tap it gently to remove any dirt from the battery charger sockets. Carefully blowing on the battery charger sockets may also help remove dirt.

- 3 Wipe the battery charger sockets with a soft dry cloth.

Cycle batteries in different sockets when recharging to wear the sockets evenly.

### If it gets wet

If the battery charger ever gets splashed with liquid, carefully shake out the liquid and dry the battery charger for 24 hours. Do not use the battery charger until it is dry.

Mains power light	What it means
 Steady green	Battery charger power is on.
Does not light up	Power adaptor is not plugged in or mains power is not available (or, if switched, is not turned on).
Battery socket light	What it means
 Steady orange	Battery module is charging.
 Steady green	Battery module is fully charged.
 Flashing orange	There is a problem with the battery module. Try a different socket. Replace the rechargeable battery module.
Does not light up	Battery module is not properly placed, is over-discharged or there is no power. Check the battery charger is powered on and reconnect the battery charger. If it still does not light up, replace the battery module.

## Using Share the Experience

See figure 12

Family members and friends can “share the experience” of bone conduction hearing. The test rod can be used by others to hear with the sound processor.

### To use the test rod

- 1 Turn on your sound processor and snap the Actuator unit onto the test rod using the tilt technique.
- 2 Place the Sound Processing unit on the ear. Hold the rod against the skull bone behind an ear. Plug both ears and listen.

## Wear an Earhook+

The Earhook+ is used when wearing the sound processor in the under-the-ear (UTE) position. It is available in different sizes and colours.

**Attach an Earhook+** See figure 13

- 1 Press firmly on the top of the Earhook with your thumb to click it off the sound processor.
- 2 Click the Earhook+ into place.
- 3 Attach and adjust the silicone band to ensure the Earhook+ sits securely.
- 4 Carefully bend the upper part into shape to suit your ear.

### Remove an Earhook+

Remove your Earhook+ by clicking it off the processor and removing the silicone band.

## Wear a Snugfit

The Snugfit holds your sound processor more securely in place than an earhook alone. It gives carers easy access to controls and accessories and is available in large, medium and tamper resistant (small).

**Attach a Snugfit** See figure 14

- 1 Press firmly on the top of the earhook with your thumb to click it off the sound processor.
- 2 Click your Snugfit into place.
- 3 Carefully bend the lower part into shape to suit your ear.

### Remove a Snugfit

- 4 Remove your Snugfit by clicking it off the processor.

**Wear a Mic Lock** See figure 15

The Mic Lock helps hold the processor securely on the ear.

### To attach the Mic Lock the first time

- 1 Push the Mic Lock band onto the processor with the tubing facing forward.
- 2 Place the processor on your ear and bring the tubing up to the earhook.
- 3 Trim the tubing to a length that will hold the processor firmly in place without discomfort. Do not cut it too short.
- 4 Soften the end of the tubing in warm water and push it onto the tip of your earhook.
- 5 Let the tubing cool before wearing.

## Attach and remove processor covers *See figure 16*

- 1 To attach, slide the cover onto the processor.
- 2 To remove, pull the cover away from the processor.

**Note:** Snugfit, Mic Lock, and Processor covers are available as accessories. For more information, ask your hearing care professional or visit [www.cochlear.com](http://www.cochlear.com)

## Replace microphone protectors *See figure 17*

Replace your microphone protectors every three months, whenever they look dirty or whenever you notice any loss in sound quality. Always replace both microphone protectors at the same time.

### Remove microphone protectors

- 1 Remove the protective cap from the Microphone Protector Removal Tool.
- 2 Push the tip of the tool into the middle of the microphone protector.
- 3 Turn the tool 90° clockwise.
- 4 Lift out the microphone protector and discard.

Repeat steps 2 to 4 to remove the other microphone protector.

### Insert new microphone protectors

*See figure 18*

- 1 Pull the Microphone Protector Applicator out of its sleeve.
- 2 With the microphone protectors facing the processor, gently push the slotted end of the applicator under the coil cable plug.
- 3 Lay the applicator over the processor, aligning the microphone protectors with the microphones.

- 4 Press the microphone protectors all the way into the microphone holes one at a time.
- 5 Remove the applicator carefully, peeling it upwards from the earhook end.

## Change the Earhook *See figure 19*

Earhooks come in large, medium and small sizes. The smallest earhook can be made tamper resistant.

- 1 Press firmly on the top of the Earhook with your thumb and click it off the sound processor.
- 2 Click the Earhook into place.

**Note:** Remove your Earhook only when necessary – it may become loose if it is removed too often.

## Assemble and change cables

*See figure 20*

Changing the cable is something you should only have to do occasionally. Using the correct method reduces the risk of accidentally damaging the cable and the cable sockets in the Sound Processing unit and the Actuator unit.

Remember to leave the Actuator unit connected when storing the sound processor in the Dry & Store unit overnight.

- 1 Hold the cable by the grips and firmly pull it straight out of the socket on the Sound Processing unit. Do not twist.
- 2 Hold the base of the plug and firmly pull the cable.
- 3 To attach a new cable, push the cable plug into the Actuator unit until it clicks into place.
- 4 Push the cable into the cable socket on the Sound Processing unit until it clicks into place.

## General processor care

Your Baha sound processor is a delicate electronic device. Follow these guidelines to keep it in proper working order:

- When not in use, turn off your sound processor and store it free from dust and dirt.
- If you will not be using your sound processor for a long period, remove the battery. Remove your sound processor before applying hair conditioners, mosquito repellent and or similar products.

## Water, sand and dirt *See figure 11*

Your sound processor is resistant to water, sand and dust. However, it is still a precision electronic device so the following precautions should be taken.

- 1 If your processor ever gets wet, dry it with a soft cloth, change the microphone protectors and place it in the dry aid kit for eight hours.
- 2 If sand or dirt ever enter the processor, shake the components carefully to remove it and change the microphone protectors.

## Warnings


### General advice

A sound processor will not restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions.

- Infrequent use of a sound processor may not enable a user to attain full benefit from it.
- The use of a sound processor is only part of hearing rehabilitation and may need to be supplemented by auditory and lip reading training.
- The patient is an intended operator of the device.
- The Model reference is Baha 5 SuperPower Sound Processor.

### Warnings

- The sound processor is a digital, electrical, medical instrument designed for specific use. As such, due care and attention must be exercised by the user at all times.
- The sound processor is not waterproof!
- Never wear it in heavy rain, in the bath or shower!
- Do not expose the sound processor to extreme temperatures. It is designed to operate within the temperature ranges +5 °C (+41 °F) to +40 °C (+104 °F). The processor should not be subjected, at any time, to temperatures below -10 °C (+14 °F) or above +55 °C (+131 °F).
- Contains small parts that may present a choking hazard.
- This product is not suitable for use in flammable and/or explosive environments.
- The sound processor and other external accessories should never be brought into a room with an MRI machine, as damage to the sound processor or the MRI equipment could occur.
- The sound processor must be removed before entering a room where an MRI scanner is located.
- If you are to undergo an MRI (Magnetic Resonance Imaging) procedure, refer to the MRI Reference Card included in the document pack.
- Portable and mobile RF (radio frequency) communications equipment can affect the performance of your sound processor.
- The sound processor is suitable for use in electromagnetic environments with mains power of typical commercial or hospital quality, and power frequency magnetic fields of typical commercial or hospital levels. It is intended for use in an electromagnetic environment in which radiated RF (radio frequency) disturbances are controlled.
- Interference may occur in the vicinity of equipment with the symbol to the right.

- Dispose of batteries and electronic items in accordance with your local regulations.
- Discard your device as electronic waste according to local regulations.
- When the wireless function is activated, the sound processor uses low-powered digitally coded transmissions in order to communicate with other wireless devices. Although unlikely, nearby electronic devices may be affected. In that case, move the sound processor away from the affected electronic device.
- When using wireless functionality and the sound processor is affected by electromagnetic interference, move away from the source of this interference.
- Be sure to deactivate wireless functionality when boarding flights.
- Turn off your wireless functionality by using the flight mode in areas where radio frequency emission is prohibited.
- Cochlear Baha wireless devices include an RF transmitter that operates in the range of 2.4 GHz–2.48 GHz.
- For wireless functionality, only use Cochlear Wireless accessories. For further guidance regarding e.g. pairing, please refer to the user guide of the relevant Cochlear Wireless accessory.
- No modification of this equipment is allowed.
- Parents and carers are advised that unsupervised use of long cables may present a risk of strangulation.
- Do not let your processor or parts entangle with any jewellery (e.g. earhook and earrings) or machinery.
- Do not attach or wear a body-worn battery configuration beneath layers of clothing in direct contact with the skin (e.g. in pocket or under a child's vest).
- Make sure all cables used by a child are securely attached to their clothing.
- Do not use a dry aid kit that has an Ultra Violet C (UVC) lamp
- The microphone protector removal tool is a sharp object. It is not recommended for use by children.
- The tamper resistant earhook fitting tool and fitting pin are sharp objects. They are not recommended for use by children
- Only use rechargeable batteries and battery chargers supplied or recommended by Cochlear. Use of other batteries or battery chargers may result in harm or injury.
- Charge rechargeable batteries before use.
- Do not allow children to replace batteries without adult supervision.
- Do not touch the battery charger contacts or allow children to use the battery charger without adult supervision.
- Do not short-circuit batteries (e.g. do not let terminals of batteries contact each other, do not place batteries loose in pockets, etc.).
- Store unused batteries in original packaging, in a clean and dry place. When processor is not in use, remove the batteries and store separately in a clean and dry place. Wipe batteries with a clean dry cloth if they become dirty.
- Battery life means the time a device will run before the rechargeable batteries must be recharged.
- Battery charge cycle is a full charge and discharge of the rechargeable battery.
- Battery lifespan means the total number of charge cycles a rechargeable battery will last before the battery life degrades to 80% of its original fullycharged capacity.
- Use only Baha 5 SuperPower Actuator units. Do not connect the sound processor to any other type of stimulating equipment such as Cochlear Implant coils.

Functions and audible/visible indicators			
Function	On the device	Audible indicator	Visible indicator
Turn ON	Connect battery	 10 beeps	 Steady light for 10 seconds
Turn OFF	Disconnect battery	None	None
Change program	Press and hold upper button	 1 beep for program 1 2 beeps for program 2 etc	 1 flash for program 1 2 flashes for program 2 etc
Change Volume	Press and release upper/lower button to increase/decrease volume	 1 beep indicates volume increase/decrease	 1 flash indicates volume increase/decrease
Low battery warning		 4 beeps 2 times	 Rapid flashes
When using Cochlear Wireless Accessory or MFi			
Change volume	With Cochlear Baha Remote Control, Cochlear Wireless Phone Clip or with an iPhone	 1 beep indicates volume increase/decrease	 1 short flash indicates volume increase/decrease
Volume limit reached		 1 long beep	 1 long flash
Enable wireless audio streaming	Press and hold lower button	 5 beeps (melody in upward scale)	 1 long flash, followed by 1, 2 or 3 short flashes
Disable wireless audio streaming	Press and hold upper button	 1 beep for program 1 2 beeps for program 2 etc	 1, 2 or 3 short flashes
Pairing confirmation		 6 beeps followed by a ripple tone (melody in upward scale)	None
Toggle between wireless accessories	Long press lower button	 5 beeps (melody in upward scale)	 Long flash, followed by one short flash
Turn ON flight mode	Insert battery and hold lower button	 10 x dual beeps	 4 x dual flashes
Mute	Mute button on Wireless Accessory or with an iPhone	 1 long beep	 1 long flash

*Note: Your hearing care professional may have disabled some or all of the audible and/or visible indicators.*

## Sound processor type designations for models included in this User Manual are:

FCC ID: QZ3BAHA5SUP, IC: 8039C-BAHA5SUP, IC model: Baha®5 SUP.

### Statement:

This device complies with Part 15 of the FCC Rules.

Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from the one in which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.
- Changes or modifications can void the user's authority to operate the equipment.

### Intended use

The Cochlear™ Baha® 5 SuperPower uses bone conduction to transmit sounds to the cochlea (inner ear). It is indicated for people with conductive hearing loss, mixed hearing loss and single sided sensorineural deafness (SSD). Furthermore it is indicated for bilateral and paediatric recipients. Fitting range up to 65 dB SNHL. It works by combining a sound processor and a small titanium implant that is placed in the skull behind the ear. The skull bone integrates with the titanium implant

through a process called osseointegration. This allows sound to be conducted via the skull bone directly to the cochlea, which improves hearing performance. The sound processor can be used together with the Baha Softband. The fitting is to be done either at a hospital, by an audiologist, or in some countries, by a hearing care professional.

### List of countries:

Not all products are available in all markets. Product availability is subject to regulatory approval in the respective markets.

The products are in compliance with the following regulatory requirements:

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**Cochlear Europe Ltd** 6 Dashwood Lang Road, Bourne Business Park, Addlestone, Surrey, KT15 2HJ, United Kingdom Tel: +44 1932 26 3400 Fax: +44 1932 26 3426

**Nihon Cochlear Co Ltd** Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo 113-0033, Japan Tel: +81 3 3817 0241 Fax: +81 3 3817 0245

**Cochlear (HK) Limited** Unit 1204, 12/F, 297-307 Hennessy Road, Wan Chai, Hong Kong Tel: +852 2530 5773 Fax: +852 2530 5183

**Cochlear Medical Device (Beijing) Co Ltd** Unit 2208-2212, Gemdale Tower B, 91 Jianguo Road, Chaoyang District, Beijing 100022, PR China Tel: +86 10 5909 7800 Fax: +86 10 5909 7900

**Cochlear Korea Ltd** 1st floor, Cheongwon Building 33, Teheran-ro 8 gil, Gangnam-gu, Seoul, Korea Tel: +82 2 533 4663 Fax: +82 2 533 8408

**Cochlear Benelux NV** Schaliënhoedreef 20 I, 2800 Mechelen, Belgium Tel: +32 1579 5511 Fax: +32 1579 5500

**Cochlear Medical Device Company India PVT Ltd** Platina Bldg, Ground Floor, Plot No. C 59, G Block, BKC, Bandra East, Mumbai 400051, India Tel: +91 22 6112 1111 Fax: +91 22 61121100

**Cochlear Italia SRL** Via Larga 33, 40138 Bologna, Italy Tel: +39 051 601 53 11 Fax: +39 051 39 20 62

**Cochlear France S.A.S.** 135 route de Saint Simon, 31035 Toulouse, France Tel: +33 5 34 63 85 85 (international) Tel: 0805 200 016 (national) Fax: +33 5 34 63 85 80

**Cochlear Nordic AB** Konstruktionsvägen 14, 435 33 Mölnlycke, Sweden Tel: +46 31 335 14 61 Fax: +46 31 335 14 60

**Cochlear Tibbi Cihazlar ve Sağlık Hizmetleri Ltd. Sti.** Cubuklu Mah. Bogazici Cad. Bogazici Plaza No: 6/1 Kavacik, 34805 Beykoz-Istanbul, Turkey Tel: +90 216 538 5900 Fax: +90 216 538 5919

**Cochlear Canada Inc** 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada Tel: +1 416 972 5082 Fax: +1 416 972 5083

**www.cochlear.com**

Cochlear Baha 5 sound processors are compatible with iPhone 6s Plus, iPhone 6s, iPhone 6 Plus, iPhone 6, iPhone SE, iPhone 5s, iPhone 5c, iPhone 5, iPad Pro, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 4, iPad mini 3, iPad mini 2, iPad mini with retina display, iPad mini and iPod touch (5th generation) running iOS 7.X or later.

"Made for iPhone" means that an electronic accessory has been designed to connect specifically to iPhone and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone may affect wireless performance.

Ardium, Baha, Baha SoftWear, BCDrive, DermaLock, EveryWear, Vistafix and Windshield are either trademarks or registered trademarks of Cochlear Bone Anchored Solutions AB. Cochlear, Nucleus, SmartSound, Snugfit, True Wireless, Hear now. And always and the elliptical logo are either trademarks or registered trademarks of Cochlear Limited. The Bluetooth word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear Limited is under license. iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. © Cochlear Bone Anchored Solutions AB 2016. All rights reserved. NOV16.

Translations of 631286-6.