93712

English Svenska Norsk


ZONE 3

MRI card

609910-2

SEPT15. Printed in Sweden.


Cochlear™ Baha® Attract System

Security control and MRI information

This person is implanted with a hearing implant system from Cochlear. The system consists of a titanium implant with a magnet implanted under the skin, and an external sound processor (1) and SP Magnet (2). The system may activate airport security metal detectors.

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Svenska

MR Safety Information

The sound processor and SP Magnet **must be removed** before entering a room where an MRI scanner is located.

Non-clinical testing has demonstrated that the BIM400 Implant Magnet, in combination with a BI300 Implant, is MR Conditional. It can be scanned safely under the following conditions. Scanning under other conditions may result in severe patient injury or device malfunction.

- Static magnetic field of 1.5 Tesla only
- Maximum spatial gradient field of 26600 Gauss/cm (266 T/m)
- Maximum switched gradient slew rate per axis of 200 mT/m/ms
- Maximum switched gradient amplitude per axis of 45 mT/m
- Maximum MR System reported whole body averaged specific absorption rate (SAR) of 2.0 W/kg (Normal Operating Mode)

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Additional instructions essential to safe use in the MR environment:

Under the scan conditions defined above, the BIM400 Implant Magnet is expected to produce a maximum temperature rise of 2.1 °C after 15 minutes of continuous scanning.

In non-clinical testing, the BIM400 Implant Magnet produced a temperature rise of less than 2.1 °C (extrapolated) at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg (extrapolated) assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla Intera, Philips Medical Systems (Software: 12.6.1.3, 2010-12-02) MR Scanner.

In non-clinical gradient-induced heating testing, the BIM400 Implant Magnet produced a temperature rise (extrapolated) of less than 4.5 °C at a time rate of change of the theoretical maximum worst-case gradient magnetic field dB/dt (extrapolated) of 200 T/s during 30 min. of continuous exposure in a test laboratory system (Pulsed Magnetic Field Generator) equivalent with a gradient system of a 1.5 Tesla MR system.

In non-clinical testing with the implant magnet in place, the image artefact caused by the device extends approximately 11.5 cm (4.5 in.) from the BIM400 Implant Magnet when imaged with a gradient echo pulse sequence and a 1.5 Tesla MRI system.

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

This person is implanted with a hearing implant system from Cochlear. The system consists of titanium implant(s) and abutment(s) penetrating the skin (1) and a sound processor (2). The sound processor can be removed before passing through a metal detector unit. If you require any additional information, contact your local Cochlear office, see www.cochlear.com for contact details.

The sound processor must be removed before entering a room where an MRI scanner is located.

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MRI Safety Information

The sound processor **must be removed** before entering a room where an MRI scanner is located.

Non-clinical testing has demonstrated that the Implant and Abutment [BI300, BIA400 & BIA300] is MR Conditional at 1.5 and 3.0 Tesla. It can be scanned safely under the following conditions.

- Static magnetic field of 1.5 Tesla and 3 Tesla only
- Maximum spatial gradient field of 3,000 Gauss/cm (30 T/m)
- Maximum MR System reported whole body averaged specific absorption rate (SAR) of 2 W/kg or maximum head averaged SAR of 3.2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the implant and abutment are expected to produce a maximum temperature rise of 1.1°C after 15 minutes of continuous scanning.

In non-clinical testing with the implant and abutment in place, the image artefact caused by the device extends approximately 1.8 cm from the implant and abutment when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. The artefact is reduced to 1.2 cm from the implant when the abutment is removed.

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